CASO CLÍNICO

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1. Stopa Z, Siewert-Gutowska M, Abed K, Szubińska-Lelonkiewicz D, Kamiński A, Fiedor P. Evaluation of the Safety and Clinical Efficacy of Allogeneic Bone Grafts in the Reconstruction of the Maxilla and Mandible. *Transplant Proc.* 2018;50:2199-2201.
2. Monje A, Chan HL, Galindo-Moreno P, Elnayef B, Suarez-Lopez del Amo F, Wang F, Wang HL. Alveolar Bone Architecture: A Systematic Review and Meta-Analysis. *J Periodontol.* 2015;86:1231-48.
3. Kuć J, Sierpińska T, Gołębiowska M. Alveolar ridge atrophy related to facial morphology in edentulous patients. *Clin Interv Aging.* 2017;12:1481-1494.
4. Hernández-Alfaro F, Sancho-Puchades M, Gujjarro-Martínez R. Total reconstruction of the atrophic maxilla with intraoral bone grafts and biomaterials: a prospective clinical study with cone beam computed tomography validation. *Int J Oral Maxillofac Implants.* 2013;28:241-51.
5. Sbordone L, Toti P, Menchini-Fabris G, Sbordone C, Guidetti F. Implant survival in maxillary and mandibular osseous onlay grafts and native bone: a 3-year clinical and computerized tomographic follow-up. *Int J Oral Maxillofac Implants.* 2009 Jul-Aug;24(4):695-703. PubMed PMID: 19885411.
6. Faot F, Marcello-Machado RM, Hermann C, Fontão FNGK. Splinted wide-short implants in the posterior region of an atrophic mandible opposed by an edentulous maxilla: immediate loading and 1-year follow-up. *Gen Dent.* 2019;67:29-33.
7. Ravidà A, Barootchi S, Askar H, Suárez-López Del Amo F, Tavelli L, Wang HL. Long-Term Effectiveness of Extra-Short (≤ 6 mm) Dental Implants: A Systematic Review. *Int J Oral Maxillofac Implants.* 2019;34:68-84.
8. Amato F. Overcoming Anatomical Limitations: The New Frontier of Implantology. *Compend Contin Educ Dent.* 2018;39:13-15.
9. Mijiritsky E, Barbu H, Lorean A, Shohat I, Danza M, Levin L. Use of Implant-Derived Minimally Invasive Sinus Floor Elevation: A Multicenter Clinical Observational Study With 12- to 65-Month Follow-Up. *J Oral Implantol.* 2016;42:343-8.
10. Franceschetti G, Trombelli L, Minenna L, Franceschetti G, Farina R. Learning Curve of a Minimally Invasive Technique for Transcrestal Sinus Floor Elevation: A Split-Group Analysis in a Prospective Case Series With Multiple Clinicians. *Implant Dent.* 2015;24:517-26.
11. Rajput N, Mohammed J. Minimally invasive transmucosal insertion and immediate provisionalization of one-piece implant in partially edentulous posterior mandible. *J Clin Diagn Res.* 2013;7:2070-3.
12. Anitua E, Alkhraisat MH. Clinical Performance of Short Dental Implants Supporting Single Crown Restoration in the Molar-Premolar Region: Cement Versus Screw Retention. *Int J Oral Maxillofac Implants.* 2019;34:969-976.
13. Anitua E, Alkhraisat MH. Fifteen-Year Follow-up of Short Dental Implants in the Completely Edentulous Jaw: Submerged Versus Nonsubmerged Healing. *Implant Dent.* 2019 Dec;28(6):551-555.
14. Anitua E, Alkhraisat MH. 15-year follow-up of short dental implants placed in the partially edentulous patient: Mandible Vs maxilla. *Ann Anat.* 2019 Mar;222:88-93.
15. Anitua E, Piñas L, Begoña L, Orive G. Long-term retrospective evaluation of short implants in the posterior areas: clinical results after 10-12 years. *J Clin Periodontol.* 2014;41:404-11.
16. Anitua E, Saracho J, Begoña L, Alkhraisat MH. Long-Term Follow-Up of 2.5-mm Narrow-Diameter Implants Supporting a Fixed Prosthesis. *Clin Implant Dent Relat Res.* 2016;18:769-77.
17. Anitua E, Errazquin JM, de Pedro J, Barrio P, Begoña L, Orive G. Clinical evaluation of Tiny® 2.5- and 3.0-mm narrow-diameter implants as definitive implants in different clinical situations: a retrospective cohort study. *Eur J Oral Implantol.* 2010;3:315-22.

18. Klein MO, Schiegnitz E, Al-Nawas B. Systematic review on success of narrow-diameter dental implants. *Int J Oral Maxillofac Implants* 2014; 29 Suppl: 43-54.
19. Ortega-Oller I, Suarez F, Galindo-Moreno P, Torrecillas-Martinez L, Monje A, Catena A, Wang HL. The influence of implant diameter on its survival: a meta-analysis based on prospective clinical trials. *J Periodontol* 2014; 85: 569-580.
20. Pommer B, Mailath-Pokorny G, Haas R, Busenlechner D, Furhauser R, Watzek G. Patients' preferences towards minimally invasive treatment alternatives for implant rehabilitation of edentulous jaws. *Eur J Oral Implantol* 2014; 7 Suppl 2: S91-109.
21. Pimentel AC, Sanches MA, Ramalho GC, Roman-Torres CV, Manzi MR, Sendyk WR. Lateralization Technique and Inferior Alveolar Nerve Transposition. *Case Rep Dent*. 2016; 2016: 4802637.