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- Al-Hashimi I, Schifter M, Lockhart P, et al. Oral lichen planus and oral lichenoid lesions: diagnostic and therapeutic considerations. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod* 2007;(Suppl. S25): e1-12.
- Mc Cartan B.E., Healy C.M. The reported prevalence of oral lichen planus: a review and critique. *J Oral Pathol Med* 2008;37:447-453.
- Scully C, El-Kom M. Lichen planus: review and update on pathogenesis. *J Oral Pathol* 1985;14:431-458.
- Roopashree MR, Gondhalekar RV, Shashikanth MC, George J, Thippesvamy SH, Shukla A. Pathogenesis of oral lichen planus- a review. *J Oral Pathol Med* 2010;39:729-34.
- George S, John SA, Anandaraj S, Issac JS, Harris A, Reshmi J. Childhood Oral Lichen Planus: Report of Two Cases. *J Dent (Tehran)*. 2015;12:374-8.
- Chandna P, Adlakha VK, Singal G, Sharma R. Pediatric oral lichen planus: review and case report. *Curr Pediatr Rev*. 2014;10:292-6.
- Pandhi D, Singal A, Bhattacharya SN. Lichen planus in childhood: a series of 316 patients. *Pediatr Dermatol*. 2014;31:59-67.
- Lavanya N, Jayanthi P, Umadevi KR, Ranganathan K. Oral lichen planus: an update on pathogenesis and treatment 2011;15:127-32.
- Sanjeev G, Chaudhari ND, Gupta A, Talanikar HV. Lichen planus- an update *Int J Pharm Biomed Sci* 2013;4:59-65.
- Carrozzo M, Thorpe R. Oral lichen planus: a review. *Minerva Stomatol* 2009;58:519-37.
- Alrashdan MS, Cirillo N, McCullough M. Oral lichen planus: a literature review and update. *Arch Dermatol Res*. 2016;308:539-51.
- Kurago ZB. Etiology and pathogenesis of oral lichen planus: an overview. *Oral Surg Oral Med Oral Pathol Oral Radiol*. 2016;122:72-80.
- Alaizari NA, Al-Maweri SA, Al-Shamiri HM, Tarakji B, Shugaa-Addin B. Hepatitis C virus infections in oral lichen planus: a systematic review and metaanalysis. *Aust Dent J*. 2016;61:282-7.
- Ma J, Zhang J, Zhang Y, Lv T, Liu J. The Magnitude of the Association between Human Papillomavirus and Oral Lichen Planus: A Meta-Analysis. *PLoS One*. 2016 29;11:e0161339.
- Sahebjamie M, Sand L, Karimi S, Biettolahi JM, Jabalameli F, Jalouli J. Prevalence of human papillomavirus in oral lichen planus in an Iranian cohort. *J Oral Maxillofac Pathol*. 2015;19:170-4.
- Scully C, Carrozzo M. Oral mucosal disease: Lichen Planus. *British Journal of oral and Maxillofacial Surgery* 2008;46:15-21.
- Mark Rodrigues M, Cherubini K, Figueiredo MA, Conclaves F. Oral lichen planus: focus no etiopathogenesis. *Archives of oral Biology* 2013;58:1057-69.
- Alberdi-Navarro J, Marichalar-Mendia X, Lartitegui-Sebastián MJ, Gainza- Cirauqui ML, Echebarria-Goikouria MA, Aguirre-Urizar JM. Histopathological characterization of the oral lichenoid disease subtypes and the relation with the clinical data. *Med Oral Patol Oral Cir Bucal*. 2017;22:e307-13.
- Carrozzo M, Thorpe R. Oral lichen planus: a review. *Minerva Stomatol* 2009;58:519-37.
- Alrashdan MS, Cirillo N, McCullough M. Oral lichen planus: a literature review and update. *Arch Dermatol Res*. 2016;308:539-51.
- Gonzalez-Moles MA, Bravo M, Gonzalez-Ruiz L, Ramos P, Gil-Montoya JA. Outcomes of oral lichen planus and oral lichenoid lesions treated with topical corticosteroid. *Oral Dis*. 2018 May;24(4):573-579.
- Expósito M, Thomsen P, Ericson L, Sennerby L, Kekholm U. Histopathologic observations on late oral implant failures. *Clinical Implant Dentistry and Related Research*. 2000;32:18-32.
- Esposito SJ, Camisa C, Morgan M. Implant retained overdentures for two patients with severe lichen planus: a clinical report. *J Prosthet Dent*. 2003;89:6-10.
- Oczakir C, Balmer S, Mericske-Stern R. Implant-prosthetic treatment for special care patients: a case series study. *Int J Prosthodont*. 2005;18:383-9.
- Reichart PA. Oral lichen planus and dental implants. Report of 3 cases. *Int J Oral Maxillofac Surg*. 2006;35:237-40.
- Hernández G, López-Pintor RM, Arriba L, Torres J, de Vicente JC. Implant treatment in patients with oral lichen planus: a prospective-controlled study. *Clin Oral Implants Res* 2012;23:726-32.
- Czerninski R, Eliezer M, Wilensky A, Soskolne A. Oral lichen planus and dental implants—a retrospective study. *Clin Implant Dent Relat Res*. 2013;15:234-42.
- Anitua E, Alkhraisat MH, Piñas L, Torre A, Eguía A. Implant-prosthetic treatment in patients with oral lichen planus: A systematic review. *Spec Care Dentist*. 2021 Jul 10.
- Anitua E, Piñas L, Escuer-Artero V, Fernández RS, Alkhraisat MH. Short dental implants in patients with oral lichen planus: a long-term follow-up. *Br J Oral Maxillofac Surg*. 2018 Apr;56(3):216-220.
- Piñas L, Alkhraisat MH, Fernández RS, Anitua E. Biological Therapy of Refractory Ulcerative Oral Lichen Planus with Plasma Rich in Growth Factors. *Am J Clin Dermatol*. 2017 Jun;18(3):429-433.
- Solakoglu Ö, Heydecke G, Amiri N, Anitua E. The use of plasma rich in growth factors (PRGF) in guided tissue regeneration and guided bone regeneration. A review of histological, immunohistochemical, histomorphometrical, radiological and clinical results in humans. *Ann Anat*. 2020 Sep;231:151528.
- Anitua E. Plasma rich in growth factors: preliminary results of use in the preparation of future sites for implants. *Int J Oral Maxillofac Implants*. 1999 Jul-Aug;14(4):529-35.
- Anitua E, Murias-Freijo A, Alkhraisat MH, Orive G. Clinical, radiographical, and histological outcomes of plasma rich in growth factors in extraction socket: a randomized controlled clinical trial. *Clin Oral Investig*. 2015 Apr;19(3):589-600.
- Yiannias JA, Rokea A, Hand JH, Pakzad AY, Rogers III RS. Relevant contact sensitivities in patients with diagnosis of oral lichen planus. *J Am Acad Dermatol* 2000;42:177-82.
- Martin MD, Broughton S, Drangsholt M. Oral Lichen Planus and Dental Materials: A Case-control Study. *Contact Dermatitis* 2003;6:331-36.
- Expósito M, Thomsen P, Ericson L, Sennerby L, Kekholm U. Histopathologic observations on late oral implant failures. *Clinical Implant Dentistry and Related Research*. 2000;32:18-32.
- Esposito SJ, Camisa C, Morgan M. Implant retained overdentures for two patients with severe lichen planus: a clinical report. *J Prosthet Dent*. 2003;89:6-10.