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1. Tan WC, Lang NP, Zwahlen M, Pjetursson BE. A systematic review of the success of sinus floor elevation and survival of implants inserted in combination with sinus floor elevation. Part II: transalveolar technique. *J Clin Periodontol*. 2008;35(8 Suppl):241-54.
2. Tavelli L, Borgonovo AE, Saleh MH, Ravidà A, Chan HL, Wang HL. Classification of Sinus Membrane Perforations Occurring During Transcresal Sinus Floor Elevation and Related Treatment. *Int J Periodontics Restorative Dent*. 2020 Jan/Feb;40(1):111-118.
3. Summers RB. A new concept in maxillary implant surgery: the osteotome technique. *Compendium*. 1994;15:152-62.
4. Bruschi GB, Crespi R, Cappare P, Gherlone E. Transcresal sinus floor elevation: a retrospective study of 46 patients up to 16 years. *Clin Implant Dent Relat Res*. 2012 Oct;14(5):759-67.
5. Chu DG, Su YC, Li YN, Zhang WW, Lyu TM, Wang LD. [Clinical applications of transcresal around detached sinus floor elevation technique using osteotomes and sequential drills with stops]. *Zhonghua Kou Qiang Yi Xue Za Zhi*. 2018 Sep 9;53(9):645-647.
6. Lundgren S, Cricchio G, Hallman M, Jungner M, Rasmussen L, Sennerby L. Sinus floor elevation procedures to enable implant placement and integration: techniques, biological aspects and clinical outcomes. *Periodontol* 2000. 2017 Feb;73(1):103-120.
7. Al-Moraissi, E. A., Altairi, N. H., Abotaleb, B., Al-Iryani, G., Helboub, E., & Alakhali, M. S. What is the most effective rehabilitation method for posterior maxillae with 4 to 8 mm of residual alveolar bone height below the maxillary sinus with implant-supported prostheses? A frequentist network meta-analysis. *Journal of Oral and Maxillofacial Surgery*, 2019;77, 70.e1-70.e33.
8. Stacchi C, Bernardello F, Spinato S, Mura R, Perelli M, Lombardi T, Troiano G, Canullo L. Intraoperative complications and early implant failure after transcresal sinus floor elevation with residual bone height ≤ 5 mm: A retrospective multicenter study. *Clin Oral Implants Res*. 2022 Aug;33(8):783-791.
9. Bernardello F, Righi D, Cosci F, Bozzoli P, Soardi CM, Spinato S. Crestal sinus lift with sequential drills and simultaneous implant placement in sites with < 5 mm of native bone: a multicenter retrospective study. *Implant Dent*. 2011 Dec;20(6):439-44.
10. Gonzalez S, Tuan MC, Ahn KM, Nowzari H. Crestal approach for maxillary sinus augmentation in patients with ≤ 4 mm of residual alveolar bone. *Clin Implant Dent Relat Res*. 2014 Dec;16(6):827-35.
11. Lafzi A, Atarabashi-Moghadam F, Amid R, Sijanivandi S. Different techniques in transalveolar maxillary sinus elevation: A literature review. *J Adv Periodontol Implant Dent*. 2021 Apr 6;13(1):35-42.
12. Anitua E, Flores J, Alkhraisat MH. Transcresal Sinus Floor Augmentation by Sequential Drilling and the Use of Plasma Rich in Growth Factors. *Int J Oral Maxillofac Implants*. 2017 May/June;32(3):e167-e173.
13. Lombardi T, Stacchi C, Berton F, Traini T, Torelli L, Di Lenarda R. Influence of Maxillary Sinus Width on New Bone Formation After Transcresal Sinus Floor Elevation: A Proof-of-Concept Prospective Cohort Study. *Implant Dent*. 2017 Apr;26(2):209-216.
14. Spinato S, Bernardello F, Galindo-Moreno P, Zaffe D. Maxillary sinus augmentation by crestal access: a retrospective study on cavity size and outcome correlation. *Clin Oral Implants Res*. 2015 Dec;26(12):1375-82.
15. Zheng X, Teng M, Zhou F, Ye J, Li G, Mo A. Influence of Maxillary Sinus Width on Transcresal Sinus Augmentation Outcomes: Radiographic Evaluation Based on Cone Beam CT. *Clin Implant Dent Relat Res*. 2016 Apr;18(2):292-300.
16. Testori T, Panda S, Clauser T, Scaini R, Zuffetti F, Capelli M, Taschieri S, Mortellaro C, Del Fabbro M. Short implants and platelet-rich fibrin for transcresal sinus floor elevation: a prospective multicenter clinical study. *J Biol Regul Homeost Agents*. 2019 Nov-Dec;33(6 Suppl. 2):121-135.
17. Anitua E, Carda C, Andia I. A novel drilling procedure and subsequent bone autograft preparation: a technical note. *Int J Oral Maxillofac Implants*. 2007 Jan-Feb;22(1):138-45. Erratum in: *Int J Oral Maxillofac Implants*. 2007 Mar-Apr;22(2):309.
18. Anitua E, Alkhraisat MH, Orive G. Novel technique for the treatment of the severely atrophied posterior mandible. *Int J Oral Maxillofac Implants*. 2013 Sep-Oct;28(5):1338-46.
19. Anitua E, Alkhraisat MH, Piñas L, Orive G. Efficacy of biologically guided implant site preparation to obtain adequate primary implant stability. *Ann Anat*. 2015 May;199:9-15.
20. Guljé FL, Raghoebar GM, Vissink A, Meijer HJA. Single crowns in the resorbed posterior maxilla supported by either 11-mm implants combined with sinus floor elevation or 6-mm implants: A 5-year randomised controlled trial. *Int J Oral Implantol (Berl)*. 2019;12(3):315-326.
21. Svezia L, Casotto F. Short Dental Implants (6 mm) Versus Standard Dental Implants (10 mm) Supporting Single Crowns in the Posterior Maxilla and/or Mandible: 2-Year Results from a Prospective Cohort Comparative Trial. *J Oral Maxillofac Res*. 2018;9(3):e4.
22. Schwartz SR. Short implants: are they a viable option in implant dentistry?. *Dent Clin North Am*. 2015;59(2):317-328.
23. Thoma DS, Haas R, Sporniak-Tutak K, Garcia A, Taylor TD, Hämmerle CHF. Randomized controlled multicentre study comparing short dental implants (6 mm) versus longer dental implants (11-15 mm) in combination with sinus floor elevation procedures: 5-Year data. *J Clin Periodontol*. 2018;45(12):1465-1474.