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1. CHANAVAZ M. Sinus grafting related to implantology. Statistical analysis of 15 years of surgical experience (1979-1994). *J Oral Implantol*. 1996; 22 (2): 119-30.
2. CHIAPASCO M, RONCHI P. Sinus lift and endosseous implants preliminary surgical and prosthetic results. *Eur J Prosthodont Restor Dent*. 1994 Sep; 3 (1): 15-21.
3. SMILER DG, JOHNSON PW, LOZADA JL, MISCH C, ROSENLICHT JL, TATUM OH JR, WAGNER JR. Sinus lift grafts and endosseous implants. Treatment of the atrophic posterior maxilla. *Dent Clin North Am*. 1992 Jan; 36 (1): 151-86; discussion 187-8.
4. WAGNER JR. A 3 1/2-year clinical evaluation of resorbable hydroxylapatite OsteoGen (HA Resorb) used for sinus lift augmentations in conjunction with the insertion of endosseous implants. *J Oral Implantol*. 1991; 17 (2): 152-64.
5. JENSEN OT. Allogeneic bone or hydroxylapatite for the sinus lift procedure? *J Oral Maxillofac Surg*. 1990 Jul; 48 (7): 771.
6. ZINNER ID, SMALL SA. Sinus-lift graft: using the maxillary sinuses to support implants. *J Am Dent Assoc*. 1996 Jan; 127 (1): 51-7.
7. SILVA LD, DE LIMA VN, FAVERANI LP, DE MENDONÇA MR, OKAMOTO R, PELLIZZER EP. Maxillary sinus lift surgery-with or without graft material? A systematic review. *Int J Oral Maxillofac Surg*. 2016 Dec; 45 (12): 1570-1576.
8. ESFAHROOD ZR, AHMADI L, KARAMI E, ASGHARI S. Short dental implants in the posterior maxilla: a review of the literature. *J Korean Assoc Oral Maxillofac Surg*. 2017 Apr; 43 (2): 70-76.
9. HADZIK J, KUBASIEWICZ ROSS P, NAWROT HADZIK I, GEDRANGE T, PITUAJ A, DOMINIAK M. Short (6 mm) and Regular Dental Implants in the Posterior Maxilla-7-Years Follow-up Study. *J Clin Med*. 2021 Mar 1; 10 (5): 940.
10. ANITUA E, ORIVE G, AGUIRRE JJ, ANDÍA I. Five-year clinical evaluation of short dental implants placed in posterior areas: a retrospective study. *J Periodontol*. 2008 Jan; 79 (1): 42-8.
11. ANITUA E, PIÑAS L, BEGOÑA L, ORIVE G. Long-term retrospective evaluation of short implants in the posterior areas: clinical results after 10-12 years. *J Clin Periodontol*. 2014 Apr; 41 (4): 404-11.
12. SUMMERS RB. A New Concept in maxillary implant surgery: the osteotome technique. *Compendium*. 1994; 15: 154-6.
13. DEL FABBRO M, CORBELLA S, WEINSTEIN T, CERESOLI V, TASCHERI S. Implant survival rates after osteotome-mediated maxillary sinus augmentation: a systematic review. *Clin Implant Dent Relat Res* 2012; 14 (Suppl 1): e159-e168.
14. ANITUA E, FLORES J, ALKHRAISAT MH. Transcrestal Sinus Floor Augmentation by Sequential Drilling and the Use of Plasma Rich in Growth Factors. *Int J Oral Maxillofac Implants*. 2017 May/June; 32 (3): e167-e173.

15. SOARDI E, COSCI F, CHECCHI V, PELLEGRINO G, BOZZOLI P, FELICE P. Radiographic analysis of a transalveolar sinus-lift technique: a multipractice retrospective study with a mean follow-up of 5 years, *Journal of Periodontology* 2013; 84: 1039–1047.
16. DEL FABBRO M, CORBELLA S, WEINSTEIN T, CERESOLI V, TASCHERI S. Implant survival rates after osteotome-mediated maxillary sinus augmentation: a systematic review *Clinical Implant Dentistry and Related Research* 2012; 14, supplement 1: e159–e168, 2012.
17. ANITUA E, PIÑAS L, ALKHRAISAT MH. Early marginal bone stability of dental implants placed in a transalveolarly augmented maxillary sinus: a controlled retrospective study of surface modification with calcium ions. *Int J Implant Dent.* 2017; 3: 49.
18. ANITUA E, ALKHRAISAT MH, PINAS L, ORIVE G. Efficacy of biologically guided implant site preparation to obtain adequate primary implant stability. *Ann Anat.* 2015; 199: 9–15.
19. ANITUA E, FLORES J, FLORES C, ALKHRAISAT MH. Long-term Outcomes of Immediate Loading of Short Implants: A Controlled Retrospective Cohort Study. *Int J Oral Maxillofac Implants.* 2016; 31: 1360-1366.
20. ANITUA E, CARDA C, ANDIA I. A novel drilling procedure and subsequent bone autograft preparation: a technical note. *Int J Oral Maxillofac Implants* 2007; 22: 138-45. Erratum in: *Int J Oral Maxillofac Implants.* 2007; 22: 309.
21. BOLLE C, FELICE P, BARAUSSE C, PISTILLI V, TRULLENQUE-ERIKSSON A, ESPOSITO M. 4 mm long vs longer implants in augmented bone in posterior atrophic jaws: 1-year post-loading results from a multicentre randomised controlled trial. *Eur J Oral Implantol.* 2018; 11 (1): 31-47.
22. GASTALDI G, FELICE P, PISTILLI V, BARAUSSE C, IPPOLITO DR, ESPOSITO M. Posterior atrophic jaws rehabilitated with prostheses supported by 5 × 5 mm implants with a nanostructured calcium-incorporated titanium surface or by longer implants in augmented bone. 3-year results from a randomised controlled trial. *Eur J Oral Implantol.* 2018; 11 (1): 49–61.
23. FAN T, LI Y, DENG WW, WU T, ZHANG W. Short Implants (5 to 8 mm) Versus Longer Implants (>8 mm) with Sinus Lifting in Atrophic Posterior Maxilla: A Meta-Analysis of RCTs. *Clin Implant Dent Relat Res.* 2017 Feb; 19 (1): 207-215.
24. GAŠPERŠIČ R, DARD M, LINDER S, OBLAK Č. The Use of 4-mm Implants Splinted to 10-mm Implants for Replacement of Multiple Missing Teeth in the Posterior Maxilla Region with Expanded Maxillary Sinus. An Observational Cases Series: Patient Characteristics and Preliminary Results. *Int J Periodontics Restorative Dent.* 2021 Mar-Apr; 41 (2): 261-268.
25. GAŠPERŠIČ R, DARD M, LINDER S, OBLAK Č. One-Year Results Assessing the Performance of Prosthetic Rehabilitations in the Posterior Maxilla Supported by 4-mm Extrashort Implants Splinted to 10-mm Implants: A Prospective Case Series. *Int J Oral Maxillofac Implants.* 2021 Mar-Apr; 36 (2): 371-378.
26. PERELLI M, ABUNDO R, CORRENTE G, SACCONI C. Short (5 and 7 mm long) porous implants in the posterior atrophic maxilla: a 5-year report of a prospective single-cohort study. *Eur J Oral Implantol.* 2012 Autumn; 5 (3): 265-72.
27. ANITUA E. Two-Step Progressive Transcrestal Sinus Augmentation Using a 4.5mm Unloaded Implant as a «Temporary Implant» in Highly Atrophic Ridge: Case Report. *Eur J Dent.* 2023 Jan 25.
28. VERRASTRO NETO A, ANDRADE R, CORRÊA MG, CASARIN RCV, CASATI MZ, PIMENTEL SP, RIBEIRO FV, CIRANO FR. The impact of different torques for the insertion of immediately loaded implants on the peri-implant levels of angiogenesis- and bone-related markers. *Int J Oral Maxillofac Surg.* 2018; 47: 651-657.
29. NORTON MR. The Influence of Low Insertion Torque on Primary Stability, Implant Survival, and Maintenance of Marginal Bone Levels: A Closed-Cohort Prospective Study. *Int J Oral Maxillofac Implants.* 2017; 32: 849-857.
30. BERARDINI M, TRISI P, SINJARI B, RUTJES AW, CAPUTI S. The Effects of High Insertion Torque Versus Low Insertion Torque on Marginal Bone Resorption and Implant Failure Rates: A Systematic Review With Meta-Analyses. *Implant Dent.* 2016; 25: 532-40.
31. WANG TM, LEE MS, WANG JS, LIN LD. The effect of implant design and bone quality on insertion torque, resonance frequency analysis, and insertion energy during implant placement in low or low- to medium-density bone. *Int J Prosthodont.* 2015; 28: 40-7.
32. ANITUA E, PIÑAS L, MURIAS A, PRADO R, TEJERO R. Effects of calcium ions on titanium surfaces for bone regeneration. *Colloids Surf B Biointerfaces.* 2015; 130: 173-81.
33. ANITUA E, PRADO R, ORIVE G, TEJERO R. Effects of calcium modified titanium implant surfaces on platelet activation, clot formation, and osseointegration. *J Biomed Mater Res A.* 2015; 103: 969-80.
34. ANITUA E, ALKHRAISAT MH, EGUÍA A. Single-crown restorations in premolar-molar regions: short (≤6.5) vs longer implants: retrospective cohort study. *Int J Implant Dent.* 2022 Oct 4; 8 (1): 40.
35. ANITUA E, FLORES J, DE RETANA SF, ALKHRAISAT MH. Single-Unit Short Implants in the Molar Region: A Retrospective Study with a Minimum 3-Year Follow-up. *Int J Periodontics Restorative Dent.* 2022 Sep-Oct; 42 (5): 683-689.
36. HERNÁNDEZ MARCOS G, HERNÁNDEZ HERRERA M, ANITUA E. Marginal Bone Loss Around Short Dental Implants Restored at Implant Level and with Transmucosal Abutment: A Retrospective Study. *Int J Oral Maxillofac Implants.* 2018 Nov/Dec; 33 (6): 1362-1367.