

JOIN THE **SQUARE** PROSTHETIC EVOLUTION

- GREATER AESTHETICS •
 - GREATER FUNCTIONALITY •
 - GREATER DURABILITY •
- in the restorations of your patients



BTI Biotechnology Institute
Tel: (+34) 945 140024 | Fax: (+34) 945 135203
pedidos@bticomercial.com
www.bti-biotechnologyinstitute.com

WE ARE bti
we are innovation


SQUARE PROSTHETIC LINE

MAXIMISES RETENTION OF
THE STRUCTURE THANKS TO ITS
4-POINT NON-ROTATIONAL BASE

PREVENTS THE BACTERIAL
FILTRATION THANKS TO ITS
HERMETIC SEAL

PROVIDES GREATER SPACE FOR THE
SOFT TISSUE AND, THEREFORE, IMPROVES
VASCULARISATION (DIRECT TO THE IMPLANT)
BIO TRANSGINGIVAL PROFILE

INCREASES THE COMPATIBILITY WITH
THE GINGIVAL FIBROBLASTS,
IMPROVING THE BIOLOGICAL SEAL
TI GOLDEN TREATMENT



**SQUARE AESTHETIC INTERFACE
OVER UNIT TRANSEPITHELIAL**

**SQUARE AESTHETIC ABUTMENT
DIRECT TO IMPLANT**